

Male Suprapubic Catheter Care Simulator

MW58


11440-000

Female Suprapubic Catheter Care Simulator

MW59

11441-000

A comprehensive trainer for suprapubic catheter changing.



Features

Changing suprapubic catheters

Trains students to gauge catheter changing skills using such metrics as catheter mobility and confirmation of urine discharge.



Catheter insertion



Confirmation of proper discharge of urine.

Skills

- Changing suprapubic catheters
- Suprapubic catheter care (securing catheter in place)
- Placement of transurethral balloon catheters.
- Enema at lateral positions


Suprapubic catheter care (securing catheter in place)

Facilitates training of securing catheter with gauze and surgical tape.



Confirmation of catheter placement

Observation of the catheter through transparent bladder wall.



Set Includes

- 1 lower torso manikin
- 1 genital organ unit
- 4 valves
(2 are incorporated in the assembled unit)
- 1 irrigation bottle
- 3 irrigation drainage tubes
- 1 instruction manual